

QUICKJETM

Beauty white.

Jet Air Polisher QUICK JET M



■ Your patients want white teeth.

The powder jet stream readily removes tobacco tar, tea tannin, stain, and plaque. With Quick Jet M, tooth whiteness and brilliance will be recovered speedily.

Compact design, handpiece type.



Easily portable spacesaving design for storage. Attachable to turbine hose and no complicated installation is required. The user-friendly designed grip fits more comfortably in your hand.



Sterilizable in an autoclave.



All components, not only the head are sterilizable in an autoclave. You can use it safely and hygienically. Inner powder passages can be cleaned from tip to toe easily and quickly.

Specifications

Available in quick-disconnect couplings of major brands.



Quick Jet M is available various models in quick-disconnect couplings of major handpiece manufacturers as well as in 2-Hole and 4-Hole connectors to meet demands of users.

More cleaning efficiency with newly designed head.

The smaller nozzle outlet enables pinpointed and efficient cleaning in a short time. The rotary head easily cleans areas hard-to-reach with conventional jet polishers.



Available in 2 colors.



Mach and Phatelus® are registered trademarks of NAKANISHI INC.

Specifications are subject to change without prior notice for improvements.

Quick Jet M comes in 2 colors.

ivory whiteblack

You can choose the preferred and suitable body color for your dental office.

Model	QUICK JET M - K	QUICK JET M - W	QUICK JET M - S	QUICK JET M - B	QUICK JET M - N	QUICK JET M - 4H	QUICK JET M - 2H
Connection	MULTIflex [®] coupling	Roto Quick [®] coupling	Sirona [®] Quick coupling	Unifix [®] coupling	Mach [®] / Phatelus [®] coupling	ISO standard 4-Hole	ISO standard 2-Hole
Powder consumption volume	4 g/min at coupling air pressure 0.4 MPa						
Air pressure	0.25 to 0.4 MPa (36 to 58 psi, 2.5 to 4.0 bar, 2.5 to 4.1 kgf/cm ²)						
Length	201 mm (7.9 inch)	200 mm (7.9 inch)	196 mm (7.7 inch)	196 mm (7.7 inch)	208 mm (8.2 inch)	207 mm (8.1 inch)	204 mm (8.0 inch)
Weight	142 g (5.0 oz.)	136 g (4.8 oz.)	137 g (4.8 oz.)	133 g (4.7 oz.)	144 g (5.1 oz.)	143 g (5.0 oz.)	143 g (5.0 oz.)
MULTIflex [®] is a registered trademark of Ka	Unifi	Unifix [®] is a registered trademark of Bien-Air Holding SA.					

MULTIflex[®] is a registered trademark of KaVo Dental GmbH. Roto Quick[®] is a registered trademark of W&H Deutschland GmbH. Sirona[®] is a registered trademark of Sirona Dental Systems GmbH.

Manufacturer : MICRON CORPORATION

2-17-7, Ikegami, Ota-ku, Tokyo, 146-0082, Japan Tel : +81 3 3755 0396, Fax : +81 3 5747 5396 http://www.micdent.com/ E-mail : contact@micdent.com

European Authorized Representative : **Advena Ltd.** Tower Business Centre, 2nd Flr., Tower Street, Swatar, BKR 4013 Malta