

## Air Scaler

# TA-200

S4H / S2H

A single multi-purpose handpiece for scaling and endodontic procedures.

### More efficient and rapid



A versatile handpiece plays double role as air scaler and endodontic instrument by changing attachments. With the vibration by compressed air, TA-200 enables much more efficient and rapid root canal cleaning and enlargement as well as tartar removal than manual treatment.

### Controllable power output



TA-200 is equipped with a power controller. The best power output can be obtained according to the purpose of use, and to the attached instrument.

Scaling

Endodontic



## Specifications

Model	TA-200-S4H	TA-200-S2H
Backend configuration	MIDWEST 4 Hole	BORDEN 2 or 3 Hole
Frequency	about 6,000 Hz	
Coupling pressure	0.3 MPa (45 psi, 3.0 bar, 3.0 kgf / cm <sup>2</sup> )	
Weight (without scaling tip)	54 g (1.9 oz)	53 g (1.9 oz)
Length (without scaling tip)	96 mm (3.8 inch)	90 mm (3.6 inch)

Specifications are subject to change without prior notice for improvements.

### Manufacturer : MICRON CORPORATION

2-17-7, Ikegami, Ota-ku, Tokyo, 146-0082, Japan  
 Tel : +81 3 3755 0396, Fax : +81 3 5747 5396  
<http://www.micdent.com/english/> E-mail : office@micdent.com

### European Authorized Representative : EMERGO EUROPE

Prinsessegracht 20 2514 AP The Hague The Netherlands  
 Tel : +31 70 345 8570, Fax : +31 70 346 7299